

# GO Pump<sup>®</sup> Rapid Recovery System



Ambulatory pain pump and catheter set  
for surgical site infusion

# The GOPump<sup>®</sup> Advantage

You take care of the patient, we'll take care of the rest

As the market leader in regional anesthesia<sup>1</sup> and provider of outsourced pharmacy compounding, B. Braun can help streamline your surgical pain pump program. Whether it's a pre-fill option through CAPS<sup>®</sup> or a 3rd party billing program through Symbios<sup>®</sup>, let us take care of the logistics and support so you can focus on what matters the most.

## B. Braun Support Services

### CAPS<sup>®</sup> GOPump Filling Service

- 22 pharmacy locations licensed in all 50 states, CAPS is registered with the FDA.
- Outsourcing allows you to free-up hospital pharmacy resources and help you drive USP <797> compliance.
- Local, reliable delivery—next day, same-day, 365 days a year.
- Ropivacaine and Bupivacaine available with 30-day shelf life.

### B. Braun Support Services

- Two dedicated sales teams covering Hospital Care and Outpatient facilities.
- Clinical and Technical telephone support team staffed by certified nurses at **800-854-6851**.
- Aesculap Academy. Provides local workshops tailored to a specific facility or group.



## Symbios Support Services

### 3rd Party Billing Service

Symbios Medical operates a pump 3rd party billing system for outpatient facilities which may help:

- Increase facility profitability.
- Reduce working capital.

### Patient Hotline Service

Symbios Medical operates a patient hotline 24/7 for troubleshooting at **866-844-5443**.



Continuous infusion of local anesthetic at the surgical site.



## Pump Simplicity

Easy to use. Easy to carry.

- Fills easily with minimal force. No need for additional equipment to make filling easier.
- No programming or rate changing necessary, just fill and go.
- Minimal patient training required.
- Convenient carrying pouch dedicated to pump size. No "one size fits all" philosophy.

## Pump Safety

Preset parameters intended to provide safety for in-patient and home use.

- Pump components do not contain natural rubber latex or DEHP.
- Durable outer shell for increased stability in home care setting. Shell can also be used as reference point to show balloon is deflating.
- Tamper evident—patient can't change anything on the pump.
- Air-eliminating filter helps prevent air from entering line.



## Items Included in Kit with GOPump®

### Choice of Catheters

Used in 100,000's of surgeries<sup>2</sup>



- Fenestrated catheter for increased anesthetic distribution. Available in 2.5", 5" and 10".



- Standard 3-port catheter.


### BIOPATCH®

- Clinical trial results showed that the use of BIOPATCH resulted in a statistically significant reduction in the incidence of CRBSI.
- BIOPATCH® Protective Disk with CHG is a hydrophilic polyurethane absorptive foam with Chlorhexidine Gluconate (CHG). The foam material absorbs up to eight times its own weight in fluid, while the CHG incorporated into the dressing inhibits bacterial growth under the dressing.<sup>3</sup>
- Chlorhexidine Gluconate is a well-known antiseptic agent with broad-spectrum antimicrobial and antifungal activity.<sup>3</sup>

For external use only. Do not use on patients with a known sensitivity to chlorhexidine gluconate.  
Refer to product labeling for complete instructions for use.

### Introducer Cannula

- Non T-Peel introducer cannula is designed to help prevent catheter site leakage and T-Peel being left behind in the patient.



## GOPump® Ordering Information

Reorder #	Description of Product(s)	Volume	Flow Rate	Catheter Type	Number of Catheters	Introducer Needle	Case Quantity
510556-BP	SP, 2.5" Cath, 150ml, 2 ml/hr, 2.75 " needle	150 ml	2ml/hr	2.5" fenestration	1	2.75"	5
510110-BP	SP, 2.5" Cath, 150ml, 2 ml/hr, 6" needle	150 ml	2ml/hr	2.5" fenestration	1	6"	5
510350-BP	SP, Epi Cath, 150ml, 2ml/hr, 2.75 " needle	150 ml	2ml/hr	Standard Catheter 3 lat eyes	1	2.75"	5
510042-BP	SP, Epid Cath, 150ml, 2ml/hr, 6" needle	150 ml	2ml/hr	Standard Catheter 3 lat eyes	1	6"	5
510112-BP	DP, 2.5" Cath, 300ml, 2 ml/hr, 6" needle	300 ml	2ml/hr	2.5" fenestration	2	6"	5
510545-BP	DP, 5" Cath, 300ml, 2 ml/hr, 6" needle	300 ml	2ml/hr	5" fenestration	2	6"	5
510201-BP	SP, 5" Cath, 300ml, 2 ml/hr, 6" needle	300 ml	2ml/hr	5" fenestration	1	6"	5
510072-BP	DP, Epid Cath, 300ml, 2ml/hr, 6" needle	300 ml	2ml/hr	Standard Catheter 3 lat eyes	2	6"	5
510401-BP	DP, 7.5" Cath, 300ml, 2ml/hr, 10" needle	300 ml	2ml/hr	7.5" fenestration	2	10"	5
510404-BP	DP, 7.5" and 5" Catheters, 300ml, 2ml/hr, 6" and 10" needle	300 ml	2ml/hr	7.5" and 5" fenestration	2	6" and 10"	5
510500-BP	DP, Epid Cath, 300ml, 2ml/hr, 2.75 " needle	300 ml	2ml/hr	Standard Catheter 3 lat eyes	2	2.75"	5
510551-BP	SP, 2.5" Cath, 300ml, 4ml/hr, 2.75" needle	300 ml	4ml/hr	2.5" fenestration	1	2.75"	5
510141-BP	DP, 10" Cath, 300ml, 2 ml/hr, 10" needle	300 ml	2ml/hr	10" fenestration	2	10"	5
510204-BP	SP, 2.5" Cath, 300ml, 4 ml/hr, 6" needle	300 ml	4ml/hr	2.5" fenestration	1	6"	5
510205-BP	SP, 5" Cath, 300ml, 4 ml/hr, 6" needle	300 ml	4ml/hr	5" fenestration	1	6"	5
510558-BP	DP, 10" Cath, 600ml, 2 ml/hr, 6" needle	600 ml	2ml/hr	10" fenestration	2	6"	5
510462-BP	DP, 5" Cath, 600ml, 2ml/hr, 6" needle	600 ml	2ml/hr	5" fenestration	2	6"	5

### Accessories

Reorder #	Description of Product(s)	Case Quantity
AN-E000	2.75" Introducer needle accessory pack	5
AN-F060	6" Introducer needle accessory pack	5
AN-F100	10" Introducer needle accessory pack	5
AP-E000	Standard Epid Cath accessory pack	5
AP-E000-275	Standard Epid Cath accessory pack w/2.75" introducer needle	5
AP-F025	2.5" fenestrated catheter accessory pack	5
AP-F025-275	2.5" fenestrated catheter accessory pack w/2.75" introducer needle	5
AP-F050	5" fenestrated catheter accessory pack	5
AP-F075	7.5" fenestrated catheter accessory pack	5
AP-F100	10" fenestrated catheter accessory pack	5

#### References:

- 1) IMS data, Q1, Q2 2010
- 2) OEM sales data on file
- 3) BIOPATCH labeling. Ethicon inc.

**Schedule a demo today.**

Call your B. Braun Representative  
at **1-800-227-2862**



**B | BRAUN**  
SHARING EXPERTISE

[www.bbraunusa.com](http://www.bbraunusa.com)